

Deluxe Decontamination Kit

Part Number: 4525-562



Introduction

The Ludlum Deluxe Decontamination Kit contains the essential equipment, labeling materials, protective clothing, decontamination solvents and cleaning supplies, needed to deal with a basic radioactive spill or routine decontamination problem that may be encountered in a nuclear medicine department or radiologic research laboratory. The 30 gallon drum serves both as a container for the kit components as well as a waste transfer and temporary storage (decay) vessel.

Specifications

Actual size & weight may vary due to options (see below) or customization of the elements of the kit.

SHIPPING SIZE (approximate): 55.9 x 55.9 x 91.4 cm
(22 x 22 x 36 in.) (L x W x H)

SHIPPING WEIGHT (approximate): 19.7 kg (43.5 lb)

Standard Kit Components

- 1 ea. 30 gallon drum (hazmat rated)
- 2 ea. disposable coveralls
- 2 pairs shoe covers
- 2 ea. N-95 mask
- 4 pairs disposable gloves
- 1 DeconGel (1 liter)
- 1 box towelettes (50/box)
- 10 ea. 30 gallon trash bags
- 1 ea. crucible tongs (18 in.)
- 2 ea. mop heads
- 1 ea. mop handle
- 1 plastic pail (3 gal.)
- 1 steel pail with lid (5 gal.)
- 1 scrub brush (6 in.)
- 1 roll (3 in. x 1000 ft) "Caution Rad Area" barricade tape
- 2 ea. laminated (7 in. x 10 in.) "Radioactive Material" signs
- 3 ea. vinyl (5 in. x 7 in.) "Radioactive Material" labels
- 1 roll absorbent paper towels